

Tendercare Tendercare Ltd, PO Box 3091
 Littlehampton, BN16 2WF
 Tel 01903 726161,
 Fax 01903 734083
 Email : info@tendercareltd.com
 www.tendercareltd.com

HARNESING ASSESSMENT FORM

Discounted prices apply when you order 4 items or more of any combination of harness/size/style.

HARNESING FEATURES

- Water & dirt-resistant
- Covered clasps
- In accordance with the CE 93/42/EEC Standard
- Non-slip grip on the interior
- Includes fastening materials
- Trim strip on adjustable strap making correct attachment even simpler

Client Name: _____

Contact Name: _____


Address: _____

Telephone No: _____




Fax No: _____

E-mail: _____

CHEST POSITIONING HARNESS'S - WATERPROOF NEOPRENE FLEXIBLE CHEST POSITIONING HARNESS



Trunk length measured, when seated, from the jugular/suprasternal notch to the top of the pelvis (iliac crest)



Shoulder width—Sizing chart (mm)	
Extra Small	230-280
Small	260 -320
Medium	320-410
Large	410-490

Discounted prices will apply when you order 4 items or more of any combination of harness/size/style.

Size (mm) & Part code	Style	A	B	C	D	E	Tick	Qty
Paed. Slim Fit KB.CHEST.TDC	Slim Fit	210	110	115	55	580		
Extra Small KB 5.71012.0	Standard	210	100	140	110	600		
Small KB 5.71012.1	Standard	300	110	210	145	600		
Small KBZ5.71022.1	With Zip	300	100	195	150	600		
Small KB 5.71014.1	Narrow/Female	300	180	220	70	600		
Medium KB 5.71012.2	Standard	350	130	230	180	600		
Medium KBZ5.71022.2	With Zip	350	130	230	180	600		
Medium KB 5.71014.2	Narrow/Female	360	240	250	70	600		
Large KB 5.71012.3	Standard	400	135	250	200	600		

FLEX CONTROL NON-STRETCH CHEST POSITIONING HARNESS


The new flex-control lets you determine the level of flexibility to suit your needs. The trim strip on the adjustable strap makes attachment even simpler.

Shoulder width—Sizing chart (mm)	
Small	260 -320
Medium	320-410
Large	410-490

Size (mm) & Part code	Style	A	B	C	D	E	Tick	Qty
Small KB 5.71032.1	Standard	300	110	210	145	600		
Small KB 5.71034.1	Narrow/Female	300	180	220	70	600		
Medium KB 5.71032.2	Standard	350	130	230	180	600		
Medium KB 5.71034.2	Narrow/Female	360	240	250	70	600		
Large KB 5.71032.3	Standard	400	135	250	200	600		

Discounted prices apply when you order 4 items or more of
any combination of harness/size/style.

H STYLE POSITIONING HARNESS - FLEX CONTROL NON-STRETCH FABRIC





Shoulder width—Sizing chart (mm)	
Small	260 -320
Medium	320-410
Large	410-490
Extra Large	490 <

Size (mm) & Part code	Style	A	B	C	D	Tick	Qty
Small KB 5.81001.1	Standard	320	50	80	80		
Medium KB 5.81001.2	Standard	400	50	100	100		
Large KB 5.81001.3	Standard	480	50	120	120		
Extra Large KB 5.81001.4	Standard	540	70	150	140		

Discounted prices apply when you order 4 items or more of
any combination of harness/size/style.

HIP POSITIONING BELTS—2 POINT HIP BELTS - SINGLE PULL ADJUSTABLE



Hip Width (HW) measured in a straight line between the top of each hip bone (iliac crest)

Measuring guide for pelvic positioning belt

Hip width—Sizing chart (mm)

Extra Small	130 - 220
Small	180 - 270
Medium	250 - 370
Large	330 - 460
Extra Large	440 <

Our pelvic belts come with two padded sections to ensure user comfort and correct pressure distribution. The nylon straps have been worked neatly into the padded sections. The padded parts are made of high-quality materials and are dirt and water-resistant.

The sizing is marked clearly on the inside of the belt and the trim strip on the nylon straps ensures correct symmetric settings. This also means you can use this to determine and check your adjustments.

SIDE RELEASE BUCKLE

Size (mm) & Part code	Style	A	B	C	D	Tick	Qty
Extra Small KB 6.70002.0	Single pull adjustable	50	175	40	600		
Small KB 6.70002.1	Single pull adjustable	70	175	50	650		
Medium KB 6.70002.2	Single pull adjustable	75	220	55	750		
Large KB 6.70002.3	Single pull adjustable	80	280	55	800		
Extra Large KB 6.70002.4	Single pull adjustable	95	330	60	800		

METAL PUSH BUTTON

Size (mm) & Part code	Style	A	B	C	D	Tick	Qty
Extra Small KB 6.70003.0	Single pull adjustable	50	175	40	600		
Small KB 6.70003.1	Single pull adjustable	70	175	50	650		
Medium KB 6.70003.2	Single pull adjustable	75	220	55	750		
Large KB 6.70003.3	Single pull adjustable	80	280	55	800		
Extra Large KB 6.70003.4	Single pull adjustable	95	330	60	800		

Discounted prices apply when you order 4 items or more of any combination of harness/size/style.

HIP POSITIONING BELTS—2 POINT HIP BELTS - SINGLE PULL ADJUSTABLE

AVIATION LATCH BUCKLE							
Size (mm) & Part code	Style	A	B	C	D	Tick	Qty
Small KB 6.70006.1	Single pull adjustable	70	175	50	650		
Medium KB 6.70006.2	Single pull adjustable	75	220	55	750		
Large KB 6.70006.3	Single pull adjustable	80	280	55	800		
Extra Large KB 6.70006.4	Single pull adjustable	95	330	60	800		

Discounted prices apply when you order 4 items or more of
any combination of harness/size/style.

HIP POSITIONING BELTS - 2 POINT HIP BELTS - DUAL PULL ADJUSTABLE


METAL PUSH BUTTON							
Size (mm) & Part code	Style	A	B	C	D	Tick	Qty
Extra Small KB 6.70004.0	Dual pull adjustable	50	175	40	600		
Small KB 6.70004.1	Dual pull adjustable	70	175	50	650		
Medium KB 6.70004.2	Dual pull adjustable	75	220	55	750		
Large KB 6.70004.3	Dual pull adjustable	80	280	55	800		
Extra Large KB 6.70004.4	Dual pull adjustable	95	330	60	800		

Discounted prices apply when you order 4 items or more of
any combination of harness/size/style.

REDUCED ACCESS PIN RELEASE							
Size (mm) & Part code	Style	A	B	C	D	Tick	Qty
Extra Small KB 6.70005.0	Dual pull adjustable	50	175	40	600		
Small KB 6.70005.1	Dual pull adjustable	70	175	50	650		
Medium KB 6.70005.2	Dual pull adjustable	75	220	55	750		
Large KB 6.70005.3	Dual pull adjustable	80	280	55	800		
Extra Large KB 6.70005.4	Dual pull adjustable	95	330	60	800		

Discounted prices apply when you order 4 items or more of
any combination of harness/size/style.

HIP POSITIONING BELTS – 4 POINT HIP BELTS - SINGLE PULL ADJUSTABLE




Hip Width (HW) measured in a straight line between the top of each hip bone (iliac crest)

Measuring guide for pelvic positioning belt

Hip width—Sizing chart (mm)

Extra Small	130 - 220
Small	180 - 270
Medium	250 - 370
Large	330 - 460
Extra Large	440 <



SIDE RELEASE BUCKLE

Size (mm) & Part code	Style	A	B	C	D	E	Tick	Qty
Extra Small KB 6.70042.0	Single pull adjustable	50	175	40	600	500		
Small KB 6.70042.1	Single pull adjustable	70	175	50	650	500		
Medium KB 6.70042.2	Single pull adjustable	75	220	55	750	600		
Large KB 6.70042.3	Single pull adjustable	80	280	55	800	600		
Extra Large KB 6.70042.4	Single pull adjustable	95	330	60	800	700		

Discounted prices apply when you order 4 items or more of any combination of harness/size/style.

HIP POSITIONING BELTS – 4 POINT HIP BELTS - DUAL PULL ADJUSTABLE

METAL PUSH BUTTON

Size (mm) & Part code	Style	A	B	C	D	E	Tick	Qty
Extra Small KB 6.70043.0	Dual pull adjustable	50	175	40	600	500		
Small KB 6.70043.1	Dual pull adjustable	70	175	50	650	500		
Medium KB 6.70043.2	Dual pull adjustable	75	220	55	750	600		
Large KB 6.70043.3	Dual pull adjustable	80	280	55	800	600		
Extra Large KB 6.70043.4	Dual pull adjustable	95	330	60	800	700		

Discounted prices apply when you order 4 items or more of any combination of harness/size/style.

HIP POSITIONING BELTS – 4 POINT HIP BELTS - DUAL PULL ADJUSTABLE

SIDE RELEASE BUCKLE

Size (mm) & Part code	Style	A	B	C	D	E	Tick	Qty
Extra Small KB 6.70044.0	Dual pull adjustable	50	175	40	600	500		
Small KB 6.70044.1	Dual pull adjustable	70	175	50	650	500		
Medium KB 6.70044.2	Dual pull adjustable	75	220	55	750	600		
Large KB 6.70044.3	Dual pull adjustable	80	280	55	800	600		


Discounted prices apply when you order 4 items or more of any combination of harness/size/style.

REDUCED ACCESS PIN RELEASE

Size (mm) & Part code	Style	A	B	C	D	E	Tick	Qty
Extra Small KB 6.70045.0	Dual pull adjustable	50	175	40	600	500		
Small KB 6.70045.1	Dual pull adjustable	70	175	50	650	500		
Medium KB 6.70045.2	Dual pull adjustable	75	220	55	750	600		
Large KB 6.70045.3	Dual pull adjustable	80	280	55	800	600		

Discounted prices apply when you order 4 items or more of any combination of harness/size/style.


ABDUCTION POSITIONING



Hip width—Sizing chart (mm)

Extra Small	130- 220
Small	180 - 270
Medium	250- 370
Large	330 - 460

Hip Width (HW) measured in straight line between the top of each hip bone (iliac crest)




ABDUCTION POSITIONING


Size (mm) & Part code	A	B	C	D	Tick	Qty
Extra Small KB 9.70001.0	190	150	70	140		
Small KB 9.70001.1	215	180	90	190		
Medium KB 9.70001.2	240	210	110	240		
Large KB 9.70001.3	270	240	130	300		

Discounted prices apply when you order 4 items or more of any combination of harness/size/style.

3 POINT PELVIC POSITIONING BELT



Hip Width (HW) measured in a straight line between the top of each hip bone (iliac crest)



Measuring guide for pelvic positioning belt



Hip width—Sizing chart (mm)	
Small	180 - 270
Medium	250 - 370
Large	330 - 460
Extra Large	440 <

3 POINT PELVIC POSITIONING BELT

Size (mm) & Part code	Style	A	B	C	D	E	F	Tick	Qty
Small KB 7.70001.1	Non Stretch	180	180	45	45	550	400		
Medium KB 7.70001.2	Non Stretch	220	220	50	50	600	450		
Large KB 7.70001.3	Non Stretch	300	300	60	60	650	500		
Extra Large KB7.70001.4	Non Stretch	400	400	80	80	750	550		
Small KB 7.70002.1	Stretch	180	180	45	45	550	400		
Medium KB 7.70002.2	Stretch	220	220	50	50	600	450		
Large KB 7.70002.3	Stretch	300	300	60	60	650	500		
Extra Large KB7.70002.4	Stretch	400	400	80	80	750	550		

Discounted prices apply when you order 4 items or more of any combination of harness/size/style.

ANKLE POSITIONING CUFFS



ANKLE CUFFS	
Ankle Circumference (A) mm	
Extra Small	< 180
Small	180-220
Medium	220-260
Large	> 260

Circumference measured over the rounded protuberance of the ankle (malleolus)

Ankle Positioning Cuffs keep the legs at the correct angle to provide a strong sitting position. If you also want to inhibit spasms, ankle positioning is an essential component in the positioning range.

Our ankle positioning device is ergonomic and provides the right amount of pressure and comfort levels. The adjustable straps mean the feet and legs still have the correct range of motion.

ANKLE POSITIONING CUFFS (PAIR)										
Size (mm) & Part code	Style	A	B	C	D	E	F	G	Tick	Qty
Extra Small KB 8.70001.0	Hook & Loop Fastening	80	65	55	45	80	210	210		
Extra Small KB 8.70002.0	Buckle Fastening	80	65	55	45	80	210	210		
Small KB 8.70001.1	Hook & Loop Fastening	110	75	60	55	110	210	210		
Small KB 8.70002.1	Buckle Fastening	110	75	60	55	110	210	210		
Medium KB 8.70001.2	Hook & Loop Fastening	140	95	65	60	140	260	260		
Medium KB 8.70002.2	Buckle Fastening	140	95	65	60	140	260	260		
Large KB 8.70001.3	Hook & Loop Fastening	170	110	65	60	170	260	260		
Large KB 8.70002.3	Buckle Fastening	170	110	65	60	170	260	260		

Discounted prices apply when you order 4 items or more of

any combination of harness/size/style.

Any Additional Information