


Adaptable to the majority of seating systems

The adaptable system

The Tilt & Fold Gas Strut Wheelbase

Tendercare

This unique wheelbase has been designed specifically for use with a wide range of special seating systems and interfaces.



	Wheelchair OPEN	Wheelchair CLOSED
Dim A	935 mm	845 mm
Dim B	690 (630) mm	360 mm
Dim C	1030 mm	670 mm





Wheelbase

One of the many features of the wheelbase is its ability to be folded compactly, thereby taking up minimal space when carried in a car, or stored in the house.

The Dual Gas Strut operated Tilt in Space facility allows it to be easily adjusted to a suitable angle for the occupant's needs. The push handles are height adjustable and reversible for the attendant's comfort and ease of use.

Interface



Moulded Seat Insert



Modular Seating System

Available in two sizes with various seat widths to help accommodate the majority of seating systems.

There is a comprehensive range of accessories available to enhance the lifestyle of individual users and their helpers.

The wheelbase has been successfully impact tested, used in conjunction with our Modular Seating System and Moulded Seat Insert.

The Wheelchair Interface System



The Seating Interface has been designed to provide a simple, quick and secure method of attaching special seating systems to wheelchairs




MODULAR SYSTEM
Tendercare



MOULDED SEAT
Tendercare



ALTERNATIVE SEATING SYSTEMS



This unique design offers a system that is easily adjustable for length and width, which also comes complete with an integral footrest location point. Therefore it is ideally suited for use with a wide range of seating systems, wheelbases and wheelchairs.

The interface is supplied with standard hooks that position the interface level with the wheelchair seat rails. Lower hooks are also available as these allow the interface to be positioned lower in the wheelbase.

The Interface is manufactured from chrome-plated steel, giving a robust long lasting construction, whilst enhancing the appearance of the wheelchair/special seat assembly and minimising maintenance requirements other than to clean it as and when necessary.

The Seating Interface has been successfully impact tested when used in conjunction with our wheelbase.

Tendercare

Tendercare Limited,
P.O. Box 3091,
Littlehampton BN16 2WF

Telephone: (01903) 726161
Facsimile: (01903) 734083
E-mail: info@tendercareltd.com
Website: www.tendercareltd.com

