



Tendercare Ltd
Unit 10, Minster Court,
Courtwick Lane,
Littlehampton, West Sussex .
BN17 7RN
Tel 01903 726161

Email: info@tendercareltd.com
Website: www.tendercareltd.com

Client Name _____

Contact Name _____

Address _____

Telephone No _____

E-mail _____

SWIFTY PUSHCHAIR

Note: The Swifty Pushchair has been impact tested and is suitable for use as a vehicle seat. It is supplied with the Tie Down Transportation Brackets fitted. The crash test was performed using the 5-point harness. This harness must be used under the belt of the passenger restraint system, when used in a wheelchair adapted vehicle. To avoid breakage always use the tilt assist when tilting chair, do not just use the push handles.

Swifty (Size 3) — Complete ☐ TH680901

Upholstery Colour Black

Maximum Carrying Weight of Occupant and all accessories 82kg
Suitable for user of height between 150-175cm

SWIFTY (Size 3) PUSHCHAIR



Complete Includes:-

Foldable Chassis with Intergrated Seat Unit, Height Adjustable Push bar, Fitted with 5-point Harness including padded shoulder straps, Belt shield pad, Breathable Fire Retardant Upholstery, Adjustable Seat Depth, Adjustable back height, Multi positional back recline 25°, Fixed seat tilt-in-space 7.2°, Height adjustable and tiltable Footrest, Tip Assist, Puncture-proof tyres, Swiveling Front Locking wheels, Quick Release Removable Rear wheels, Foot parking brake, Set of 4 Tie Down Transportation Brackets fitted.

SPECIFICATIONS

Specifications	Swifty 3
Seat Depth	370 - 460mm
Seat Width	410mm
Back Height	680 - 830mm
Shoulder Height	490 - 650mm
Lower Leg Length	250 - 450mm
Hip Angle	90° -115°
Foot Angle	90°
Measurement footrest (w x d)	295 x 230mm
Seat Tilt	7.2°
Total (l x w x h)	1130 x 6900 x 900mm
Folded Measurements	600 x 690 x 900mm
Weight	24.5kg
Max Weight Capacity	82kg

FACTORY FITTED UPGRADES

Push Handle mounted attendant brakes - Hub* ☐ SN6967
* Standard foot operated brake is also fitted


Pneumatic Wheels Upgrade (2 x front and 2 x rear wheels) ☐ TH6966

POSTURAL ACCESSORY


Belt shield pad including magnetic lock ☐ TH434
(To help prevent unlocking of harness buckle)

LATERAL SUPPORTS


Lateral Head Supports




Lateral Trunk Supports



Lateral Trunk Supports (swing away)




(All Dimensions in mm)		✓	A	B
Lateral Head Supports (pair)	TH6955	☐	150	120
Lateral Trunk Supports (pair)	TH6956	☐	170	120
Lateral Trunk supports (pair) (Swing Away)	TH6957	☐	A1 = 150 A2 = 110 Total A=260	110

POMMEL					LEG AND FOOT POSITIONING											
 includes fixing bracket					Toe Straps- Pair (Cam Buckle) ☐ TS6591											
FRAME ACCESSORY																
<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 25%;">(All Dimensions in mm)</td> <td style="width: 5%; text-align: center;">✓</td> <td style="width: 10%;">A</td> <td style="width: 10%;">B</td> <td style="width: 10%;">C</td> </tr> <tr> <td>Pommel (Size 3) TH4121</td> <td style="text-align: center;">☐</td> <td style="text-align: center;">130</td> <td style="text-align: center;">90</td> <td style="text-align: center;">110</td> </tr> </table>					(All Dimensions in mm)	✓	A	B	C	Pommel (Size 3) TH4121	☐	130	90	110	Frame Padding ☐ TH6965	
(All Dimensions in mm)	✓	A	B	C												
Pommel (Size 3) TH4121	☐	130	90	110												

LIFESTYLE ACCESSORIES			
Sun canopy - (includes fixing brackets)	☐	TH6963	Incontinence pad (fits on top of seat base cover) ☐ TH6961 it is possible to insert other commercially available absorbent pads into the incontinence seat pad to increase the absorption capacity if required.
Rain Cape / Poncho - fits over user and pushchair frame (lateral head supports and Sun canopy do not fit under the poncho hood)	☐	TH6964	
Shopping Basket (Max carrying capacity 5kg) Basket needs to be removed for folding of chair	☐	TH6970	Under seat storage bag with zip - (Max carrying weight 5kg) ☐ TH6971 Bag needs to be removed for folding of chair

PATIENT MEASUREMENT RECORD	
Name: _____	D.O.B: _____
Diagnosis: _____	
Current Equipment: _____	
Height: _____ Weight: _____	
Date Measurements Taken: _____	
By Whom: _____	

Please state actual client measurements in mm.	
1 Top of head to seat: _____ 2 Top of shoulder to seat: _____ 3 Under armpit to seat: _____ (Axilla Height) 4 Actual sitting depth: _____ 5 Shoulder width: _____ 6 Chest Width: _____ (Arm Pit to Arm Pit) 7 Hip width: _____ 8 Left leg drop: _____ 9 Right leg drop: _____	

INFORMATION
<i>The Swifty pushchair complies with the requirements of ISO 7176-19 and is designed to be secured facing forward.</i>
ANY ADDITIONAL INFORMATION / SPECIAL REQUESTS
We make our own accessories. If you need something more bespoke/customised, please provide details below. All bespoke items will be priced depending on requirements.